

Vios™ Aerosol Delivery System - Adult Configurations

STANDARD

#310F35-LCS with PARI LC® Sprint NDC# 83490-0310-58

#310F83-LC+ with PARI LC® Plus NDC# 83490-0310-67



PARI LC® Sprint Reusable Nebulizer and backup reusable nebulizer

(OR)



PARI LC® Plus Reusable Nebulizer and backup reusable nebulizer



Timestrip® - 6 Month Reminder



Vios™ Premium Packaging



Vios™ Compressor (Blue)



Vios™ Experience Instructions

DELUXE

#310F35-LSD with PARI LC® Sprint NDC# 83490-0310-59

#310F81-LCD with PARI LC® Plus NDC# 83490-0310-64

INCLUDES ALL STANDARD COMPONENTS AND:



Dual Compartment Carrying Case



Vios™ Experience DVD Video

Questions? 1.800.FAST.NEB (327.8632)
www.PARI.com



Vios™ Compressor Technical Specifications:

Maximum Pressure of Compressor = 34 PSI (2.3 bar)

Maximum Liter Flow of Compressor = 8.5 LPM

Nominal Operating Pressure¹ = 18 PSI (1.2 bar)

Nominal Air Flow¹ = 4.5 LPM

Weight = 3 lbs (1300 gms)

Dimensions (Length, Width, Height) = 6.5 x 6.5 x 3.75 in

Warranty = 5 years, 6 years user activated

HCPCS Billing Codes:

Nebulizer with Compressor = E0570

Reusable Nebulizer = A7005

¹When used with a PARI LC® reusable nebulizer.



INTRODUCING



AEROSOL DELIVERY SYSTEM

From the leader in clinically proven aerosol delivery



PARI nebulizers are the most widely used nebulizers by pharmaceutical companies in clinical trials for new medications compared to any other brand.



WHY VIOS™?



“ Because there is a difference in aerosol delivery systems and Vios™ outperforms the others. ”



“ Because with Vios™, I am confident that my patients will receive an effective treatment dose. ”



“ Because my Vios™ treatments are fast, I can spend more time with my grandchildren. ”

The Vios™ Aerosol Delivery System is designed from PARI's long legacy of clinically proven products and consists of two important components: a Vios™ compressor and a PARI nebulizer. Vios™ was inspired by the Latin and Greek words for "Life" (Vita + Bios = Vios). Our goal is to improve quality of "Life" and patient lifestyle through efficient aerosol delivery. The Vios™ delivers consistent particle size and fast treatments resulting in targeted delivery to the lungs.

Consistent Aerosol Delivery

EXCLUSIVE TO PARI, **TIMESTRIP®** IS AN EASY TO USE VISUAL INDICATOR THAT MONITORS TIME ELAPSED FROM A PATIENT'S FIRST TREATMENT TO 6 MONTHS. NO MORE GUESS WORK OR HASSLE WITH CALENDARS - TIMESTRIP® CLEARLY INSTRUCTS PATIENTS TO REPLACE THEIR PARI REUSABLE NEBULIZER EVERY 6 MONTHS.



- Easy to use and activate
- Great for Patient Reminder programs
- Benefits of replacing PARI nebulizers:
 1. Consistent medication delivery
 2. Fast treatment times
 3. Reduces the risk for bacterial contamination